



ADJUSTABLE CRUTCHES BRACKET

INSTRUCTION MANUAL



Trying to carry a walking stick, crutches, or even an umbrella while using your wheelchair or scooter or getting around ambulatory, can get awkward fast.

This clip-on holder gives you a simple way to keep those essentials with you, without needing to balance them on your lap or juggle them in your hands. It attaches directly to your chair or scooter frame, so your walking aid or umbrella stays secure, within reach, and ready when you need it.

The two-piece design lets you position the clips at the height and spacing that works best for you and your setup. Once in place, it holds firmly, so you can head out knowing everything is where it should be.

Great for:

- Wheelchair and mobility scooter users
- Carrying walking sticks or canes
- Holding crutches securely while moving
- Ambulatory users who switch between walking and mobility aids
- Keeping an umbrella within easy reach
- Reducing the need to carry items by hand
- Everyday outings, appointments, and getting around

Features:

- Secure clip design to hold items in place
- Two-piece setup for adjustable positioning
- Clips can be placed at your preferred height and spacing
- Attaches directly to wheelchair or scooter frame
- Lightweight but sturdy construction
- Simple black design that blends in with most setups

Image Description: A black plastic clip-on holder designed to attach to a wheelchair or mobility scooter frame. It features two separate pieces: a top clip that wraps around the pole and a bottom cup-style holder that supports the base of a walking stick, crutch, or umbrella. Images show the clips from multiple angles, including how they open, attach, and hold items securely in place.



ADJUSTABLE CRUTCHES BRACKET

INSTRUCTION MANUAL

Details:

Material: ABS Plastic

Size: please see part details below...

Top piece: 7.5cm (L) x 3.5cm (W) x 3.5cm (D)

Bottom piece: 11.5cm (L) x 5.5cm (W) x 3.5cm (D)

Fits Crutches: Walking sticks or crutches approx. 1.5cm to 2cm diameter

Fits Pole: approx. 2.5cm diameter

Weight: 50 grams

Package includes: 1x Crutch Holder set (top clip + bottom holder)

Age Guidance: Suitable for all ages once installed. Adult installation recommended.

How to Use:

- Position the clips on your wheelchair or scooter frame where they suit your setup
- Secure the top clip firmly around the pole
- Attach the bottom holder at a height that supports the base of your item
- Place your walking stick, crutch, or umbrella into the holder
- Check everything feels secure before heading out
- Adjust positioning if needed for comfort or fit

Care Instructions:

- Wipe clean with a damp cloth as needed
- Avoid over-tightening to prevent damage to the clip
- Store safely if removed from your device
- Check regularly to ensure clips remain secure

Important Note:

- Fit may vary slightly depending on your wheelchair or scooter frame size. If unsure, feel free to get in touch and we can help you check compatibility.