



ADJUSTABLE TABLETOP IV POLE

INSTRUCTION MANUAL



Trying to create a setup away from home can quickly turn into cords, bags, tubes, and equipment balancing awkwardly off tables, beds, chairs, or whatever happens to be nearby. This foldable travel clamp IV pole gives you a portable place to hang and organise equipment wherever you need it, without needing permanent mounting or bulky medical stands taking over the room.

Designed with a clamp-on base, telescoping height adjustment, and collapsible lightweight frame, it is easy to pack down for travel, temporary setups, bedside use, wheelchair spaces, hotel stays, events, gaming setups, or everyday support needs at home. The adjustable pole helps bring equipment closer and more accessible while reducing the need to constantly reach, bend, balance, or improvise. Because accessibility equipment should be able to move with your life, not keep you stuck in one place.

Great for:

- Bedside feeding or hydration setups
- Travel and hotel stays
- Wheelchair tray or side table setups
- Temporary home care spaces
- Gaming and streaming setup
- Craft, sewing, or hobby organisation
- Holding lightweight medical accessories
- Campervan or portable living setups
- Folding down for transport or storage
- Clamping onto desks, tables, or shelves
- Support workers creating flexible setups
- Small spaces where full IV stands are impractical

Features:

- Lightweight portable design
- Telescoping adjustable height pole
- Foldable collapsing structure for storage
- Clamp-on mounting system



ADJUSTABLE TABLETOP IV POLE

INSTRUCTION MANUAL

- Easy-adjust height clip
- Packs down smaller for transport
- No permanent installation require
- Easy to reposition between space
- Compact black frame with blue adjustment points
- Dual hanging hooks at the top

Image Description: A portable black clamp-style IV pole with adjustable telescoping height and dual hanging hooks at the top. The pole attaches to desks, tables, trays, or flat surfaces using a screw clamp base. Several images show close-up views of the clamp mechanism, height adjustment clips, and the lightweight foldable frame designed for travel and portable setups.

Details:

Material: Stainless Steel, Aluminium Alloy, ABS Plastic

Size: Please see sizing below...

Pole Size: 40cm - 79cm (L) x 20cm (W) x 2.5cm (diameter)

Clamp Size: 7.5cm (L) x 4cm (W) x 5.5cm (D)

Full Size: 48.5cm - 87.5cm (L) x 20cm (W) x 2.5cm (diameter)

Travel Size: 56cm (L) x 8cm (W)

Weight: 1kg Approx.

Package Contains: 1x Tabletop IV Pole

Age Guidance: Intended for teens and adults. This is not a toy.

How to Use:

1. Attach the clamp securely onto a stable surface such as a desk, tray, bedside table, or shelf
2. Adjust the pole height using the telescoping clips
3. Hang lightweight compatible equipment or accessories from the top hooks
4. Reposition or move between spaces as needed
5. Collapse and fold down the pole for transport or storage

Care Instructions:

- Wipe clean with a damp cloth
- Avoid excessive weight or overloading the hooks
- Store in a dry place when not in use
- Check clamp tightness regularly before use

Important Notes:

- Surface shown in images is not included
- Ensure the clamp is attached securely before hanging equipment
- Designed for lightweight hanging support only
- Not intended for critical medical use or unsupported heavy loads