



Getting the right dose from tablets or supplements can be difficult when pills are small, uneven, or hard to split safely by hand.

This universal pill cutter is designed to help create more accurate, consistent splits for pills, vitamins, and supplements. The enclosed design and guided chambers help hold tablets in place, while the cutting mechanism allows for cleaner, more controlled division without needing knives or improvised tools.

Compact and easy to store, it is designed for everyday home use, travel, or keeping in a medication kit.

Great for:

- Splitting tablets and supplements into smaller doses
- Supporting accurate medication routines
- Reducing stress around manual pill cutting
- Arthritis, reduced grip strength, or hand coordination needs
- Elderly users managing daily medication
- Neurological or movement-related conditions support
- Travel and on-the-go medication management
- Pet medication dosing at home

Features:

- Multi-chamber pill positioning system
- Enclosed blade cutting design
- Helps reduce crushing or uneven splitting
- Lightweight and portable
- Travel-friendly size
- Safer alternative to knives or manual cutting tools

Image Description: A compact pill cutter device with a red and grey plastic body, featuring multiple circular compartments for different pill sizes and a hinged cutting lid. The images show the device open and closed, highlighting the pill alignment chambers and internal cutting blade design.

**Details:**

Material: ABS Plastic + Stainless Steel

Size: 6.8cm x 10.5cm x 4.3cm

Weight: 35 grams

Package Contains: 1x Universal Pill Cutter

Age Guidance: Recommended for adults. Use with supervision where needed for safety or medication support. Not a toy.

How to Use:

- Place tablet into the appropriate chamber
- Close lid securely to align pill
- Press down to split tablet
- Clean after use if required
- Store safely when not in use

Care Instructions:

- Wipe clean with a dry or lightly damp cloth
- Ensure fully dry before storing
- Keep blade area closed when not in use
- Store out of reach of children

Please Note:

- Designed for tablets that are suitable for splitting
- Not all medications are appropriate to cut – follow medical guidance
- Use care when handling blade mechanism