

## WALKING STICK AND CRUTCHES CLIP



### **Description:**

Keeping mobility aids together shouldn't be another thing to manage.

These crutch connector clips make it easy to keep your crutches paired together when you're not using them, so they're easier to carry, store, and grab when you need them. Using two clips, one higher up and one lower down, helps hold crutches neatly together without slipping or tangling.

They're also handy when you use more than one mobility aid. You can clip a walking stick or crutch onto a walking frame, giving you access to both options while you're out and about. It's a small, simple solution that supports independence, flexibility, and doing what works best for you.

### **Great for:**

- 📌 Keeping crutches securely paired when not in use
- 👤 Preventing one crutch being misplaced or left behind
- 🧳 Travel, storage, and transport
- 🏠 Keeping mobility aids tidy at home
- 👤🛠️ Clipping a walking stick or crutch to a walking frame
- 👤 Supporting independence and choice in how you move

### **Features:**

- 🌀 Sold as a set of two clips for better stability
- 📌 Designed to connect crutches, walking sticks, or poles
- 👤 Easy snap-on design with no tools required
- 🧳 Compact and lightweight for everyday use
- 🔄 Flexible use depending on your mobility setup

**Details:**

**Material:** ABS Plastic

**Size:** 6cm (L) x 2.2cm (W) x 1cm (D)

**Weight:** 10 grams

**Package includes:** 2x Crutches Clips

**Age Guidance:** Suitable for ages 3+ with adult supervision.

**How to Use:**

-  Clip one connector onto the upper section of a crutch or pole
-  Clip the second connector lower down for added stability
-  Press the second crutch or pole into place until secure
-  To attach to a walking frame, clip one side onto the frame and the other onto your stick or crutch
-  Adjust placement to suit what feels secure and comfortable for you

**Care Instructions:**

-  Wipe clean with a damp cloth
-  Do not force open beyond intended use
-  Store somewhere easy to reach when not in use